


Fabrum Solutions Pulse Tube Cryocooler

Model: Evercold PTC1000

Product Description V.1.0

For safe and correct use, please read the operation manual thoroughly before operating the Cryogenic refrigerator.



Manufacturer:

Fabrum Solutions Ltd

PHONE +64 3 982 2810, FAX +64 3 982 2811

ADDRESS 128a Waterloo Road, Hornby, Christchurch, New Zealand 8042

www.fabrum.com

Document Control- Version 1.0
Released: February 2016



Contents

| | |
|---|---|
| 1. Technical Specification - Brief | 3 |
| 2. Geometry and Components..... | 3 |
| 3. Condenser Layout | 5 |
| 4. Heat-lift Performance | 5 |
| 5. Air Cooled Addition | 5 |
| 6. Installation Guideline | 6 |
| 7. Operating Guidelines | 6 |
| 8. Fault Indication | 6 |
| 9. Maintenance Schedule..... | 7 |
| 9.1 Oil/filter Change..... | 7 |
| 9.2 Inspect vibration isolator | 7 |
| 9.3 Coolant pressure check..... | 7 |
| 9.4 Helium top up | 7 |
| 9.5 Replacement of Seals and wear components..... | 7 |
| 10. Compliance..... | 7 |
| 11. Warranty statement..... | 8 |

1. Technical Specification - Brief

| | |
|----------------------------|-----------------------------------|
| Model | PTC1000 |
| Swept volume | 1000 cm ³ |
| Power @ 77 K | 1300 W |
| No load temperature | 43 K |
| Motor Power @ 77 K | 25 kW |
| Weight | 1300 kg |
| Dimensions | 1.8 (L) x 1.2 (W) x 2.1 (H) |
| Noise at 3 m | 88 dBA |
| Frequency | 50/60 Hz |
| Current rating | 63 A |

| | |
|---|---------------------------------|
| Charge Pressure | 25 Bar |
| Pressure amplitude | 5 Bar |
| Ambient temperature | 0-35°C |
| Heat rejection | Water/Air |
| Cooling water flow rate (In water cooled configuration) | 40 l/min Approx. 3 Bar ΔP |
| For closed loop air cooled configuration | |
| Coolant pump power | 1.1 kW |
| Radiator fan power | 1.4 kW |
| All specifications subject to change. Cooling power are given for water cooled configuration with 15°C water temperature. Dimensions and weights are for the cryocooler only. | |

2. Geometry and Components




Figure 1 PTC1000 geometry




Figure 2 Component description

Proprietary Components

| | |
|-------------|---------------------|
| 30 kW Motor | ABB-M3AA200MLA-2-B5 |
|-------------|---------------------|

| | |
|----------------------|------------------|
| Variable speed drive | ACS580-01-072A-4 |
|----------------------|------------------|

(not shown)

3. Condenser Layout

The condenser design depicted in Figure 3 has been optimised for use with a PSA feeding gas to the inlet port. In this configuration, a working dewar is used to collect the liquid, then transferred to storage and customer dewars. The condensers are within the cryostat which is under static vacuum.

All wetted surfaces are austenitic stainless steel and copper which have good compatibility with LN2, LO2 and LNG.




Figure 3 Condenser layout

4. Heat-lift Performance




Figure 4 Cooling performance

The efficiency in Figure 4 is calculated using motor power.




Figure 5 Coolant temperature dependence.

Figure 5 Shows the Cryocoolers derating of 4.5 W at 77 K per degree rise in cooling media temperature. This data is based on good coolant flow of 40 L/min.

5. Air Cooled Addition




Figure 6 Air Cooling system.

Mass: 120 kg (dry)

Figure 6 shows the additional equipment required for air cooled operation. Separation from the cryocooler is trivial if required, alternately a frame can be supplied to support both the cryocooler, cooling system, controls and a dewar if required.

6. Installation Guideline

The unit must not be installed by the user. A qualified technician must install the unit.

- The unit must be mounted to an appropriate frame. Custom frames can be supplied to your specifications.
- This unit should be mounted on a level surface.
- For any bespoke applications please consult with the manufacturer.
- Direct sunlight (UV radiation) exposure should be avoided as it can damage some plastics used in this unit.
- Install only on robust concrete floor or similar.
- Every cryocooler unit is tested at the factory and carefully checked before shipment. However, independent to the checks made at the factory, the unit could be damaged during transport and therefore needs to be re-inspected.
- The unit must be installed indoors with adequate shelter.

7. Operating Guidelines

| Electrical connection | |
|-----------------------------|---|
| Voltage | 3-phase, $U_{N4}=380$ to 480 V, +10%/-15% |
| Frequency | 50/60 Hz $\pm 5\%$ |
| Power Factor | $\cos\phi=90$ |
| Environmental limits | |
| Ambient temperature | 0 to +35°C standard 0 to +45°C by request |
| Altitude | 0 to 1,000 m without derating 1,000 to 4,000 m with derating of 1%/100 m |
| Relative humidity | 5 to 95%, no condensation allowed. |
| Contamination levels | No conductive dust allowed |

STARTUP PROCEDURE:

- Reset the emergency stop button by twisting the knob if necessary.
- Turn switch on Cryocooler electrical cabinet to ON.

TYPICAL OPERATING LEVELS

- Oil level: 2/3 way up the sight gauge
- Cooling system pressure (outflow from pump): 65PSI
- Oil and machine temperature < 60°C.
- Helium gas pressure - 25bar.
- Helium gas pressure amplitude - 5 bar.
- Running speed 48Hz.

SHUT-DOWN PROCEDURE:

- Turn switch on liquefier electrical cabinet to OFF.

EMERGENCY SHUTDOWN PROCEDURE:

- Press emergency stop button
- Isolate power to the machine.

8. Fault Indication

| Interlock Monitor | Fault |
|-------------------|--|
| light #1 | The Oil pressure is <5 bar after 15 s of running |
| light #2 | Gas Pressure is < 5 bar |
| light #3 | Top Plate temperature > 70°C |
| light #4 | Machine temperature >70°C |
| light #5 | Water pressure (pump inlet) < 5PSI |
| light #6 | Emergency stop- pressed |

9. Maintenance Schedule

All maintenance must be performed with the unit off.

| Period | Work to be done | Chapter |
|-------------------------|--|---------|
| Annually | Oil/filter Change | 9.1 |
| Annually | Visually inspect vibration isolators | 9.2 |
| Annually | Check coolant pressure | 9.3 |
| Annually | Helium gas top up | 9.4 |
| 4 years (40, 000 hours) | Replacement of seals and wear components | 9.5 |

9.1 Oil/filter Change

Drain oil from oil reservoir by removing the plug near the base. Drain into a suitable container for disposal. Re-install plug into base of reservoir ready for refill.

Fill oil reservoir through the cap on the top of the reservoir with 35 litres of Hydraulic oil. The oil level should be to the top of the level gauge on the front face of the oil reservoir.

Remove 2x spin on oil filter from the unit and replace.

| Required equipment | |
|--------------------|-----------------|
| Hydraulic Oil | Shell Tellus 68 |
| Spin on oil filter | BE-10-18 filter |

9.2 Inspect vibration isolator

Visually inspect all 4 vibration isolators at the base of the machine. If the rubber is damaged these must be replaced.

9.3 Coolant pressure check

Note: This applies to closed loop air cooled unit only.

The static coolant pressure must read between 15 and 20 PSI (100 to 140 kPa). If the gauge reads outside these limits, please contact the manufacture.

9.4 Helium top up

If the running helium pressure is below 24 bar, a top up is required.

A top up can only be performed once the unit has equalized to room temperature.

Connect to Ultra High Purity Helium supply (99.999% purity) via gas spring fill valve. Ensure helium transfer valve is fully open. Fill with helium slowly to **22 bar** (approx. 1 bar/min), as the transfer line restricts helium flow between the gas spring and pulse tube. A room temperature fill to 22 bar gives 25 bar average running pressure.



9.5 Replacement of Seals and wear components

After 4 years of operation we recommend replacement of (PWG) seals, main bearings, piston seals, diaphragm, regenerator mesh (PT), and inertance tube (PT).

Safe and reliable operation of this unit is guaranteed only with Fabrum Solutions spare parts.

10. Compliance

All components within the pressure volume envelope have been designed to meet or exceed ASME VIII DIVISION 1 2013.

We are working together with a compliance consultant and Conformance UK to get CE approval for this unit. The pressure volume product is over 200 barL and therefore external certification is required. All electronics and controls are commercially available with CE marking.

11. Warranty statement

- a) Full warranty on product for four (4) years of operation, from date of commissioning, subject to noted exclusions.
- b) Warranty covers replacement parts in that four year period, and provision of a technician on site for first 12 months (labour and disbursements not covered post 12 months).
- c) Return to base freight is covered one way if return is required.
- d) Not covered for normal wear and tear as specified in the preventative maintenance schedule.
- e) Not covered for over temperature operation, tampering, or misuse.
- f) Safe and reliable operation of this unit is guaranteed only with Fabrum Solutions spare parts.
- g) Fabrum Solutions will carry spare Pulse Tubes and at least one 'ready for dispatch' replacement part for the Diaphragm Pressure Wave Generator assembly (piston seals, bearings, cranks, master and slave piston and cylinder).